

## Division of Nuclear Medicine Procedure / Protocol University Hospital and The American Center

PARATHYROID SPECT PROCESSING Tc99m  
UPDATED: JUNE 2016

CPT CODE: N/A

**Recon / Reformat Processing and Screen Captures for SPECT/CT**

Reconstruction Options:	Infinia 1 (Rm E)	Infinia 2/3 & NM640 with Xeleris 3.x (RM A, C and TAC)
Reconstruction Type	OSEM/MLEM	OSEM
Number of OSEM Iterations	OSEM 2 MLEM 0	2
Max Number of Subsets	10	10
<b>Correction Options:</b>		
Corrections		Resolutions recovery correction & Attenuation Correction
Additional Set		Resolution Recovery Correction and Attenuation Correction
<b>Filters:</b>		
Pre-Filter	Butterworth	
Pre Critical Frequency	0.50	
Pre Power	10	
Post-Filter	Hann	Butterworth
Critical Frequency	1.2	0.50
Power	10.0	10.0

Infinia 1 with Xeleris 2.x (Rm E)	Infinia 2/3 & NM640 with Xeleris 3.x (RM A, C and TAC)
1. Select the Patient and the following files	1. Select the Patient and the following files
a. TOMOearlyHWKY	a. TOMOearlyHWKY
b. CT TOMOearlyHWKY	b. CT TOMOearlyHWKY
c. ATT MAP TOMOearlyHWKY	c. ATT MAP TOMOearlyHWKY
2. Run VOLUMETRIX FOR HAWKEYE PARATHYROID processing icon	2. Run Parathyroid MI processing icon
	a. Click Original
	b. Click Proceed
	c. Click Passed
3. Select NM Transaxials	3. Select NM Transaxials
	a. Change SUM:1 and STEP:1
a. Create an 8 x 8 display grid	b. Create an 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	c. Center area of interest within the grid and adjust intensity accordingly
c. Annotate EARLY TRANS on display screen	d. Annotate EARLY TRANS on display screen
d. SCREENCAP and save as EARLY NM TRANSAXIALS	e. SCREENCAP and save as EARLY NM TRANSAXIALS

4. Repeat step 3 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately	4. Repeat step 3 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately
5. Select <b>FUSED Transaxials</b>	5. Select <b>FUSED Transaxials</b>
a. Create an 8 x 8 display grid	a. Create an 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate <b>EARLY TRANS</b> on display screen	c. Annotate <b>EARLY TRANS</b> on display screen
d. <b>SCREENCAP</b> and save as <b>EARLY FUSED TRANSAXIALS</b>	d. <b>SCREENCAP</b> and save as <b>EARLY FUSED TRANSAXIALS</b>
6. Repeat step 5 by selecting <b>FUSED CORONALS AND FUSED SAGITTALS</b> . Annotate appropriately	6. Repeat step 5 by selecting <b>FUSED CORONALS AND FUSED SAGITTALS</b> . Annotate appropriately
7. Click <b>File</b> and <b>Save MIP</b>	7. Click <b>File</b> and <b>Save MIP</b>
8. Click <b>File</b> then <b>Save &amp; Exit</b>	8. Click <b>File</b> then <b>Save &amp; Exit</b>
9. Select MIP and rename to <b>EARLY MIP</b>	9. Select MIP and rename to <b>EARLY MIP</b>
10. Select <b>CT TOMOearlyHWKY</b>	10. Select <b>CT TOMOearlyHWKY</b>
11. Run <b>Convert CT to Hounsfield Units</b> processing icon	11. Run <b>Convert CT to Hounsfield Units</b> processing icon
12. Click <b>File</b> and <b>Exit</b>	12. Click <b>File</b> and <b>Exit</b>
13. Select <b>Tomolate</b>	13. Select <b>Tomolate</b>
14. Run <b>VOLUMETRIX FOR HAWKEYE PARATHYROID</b> processing icon	14. Run <b>VOLUMETRIX FOR HAWKEYE PARATHYROID</b> processing icon
15. Click <b>Resume</b>	15. Click <b>Resume</b>
16. Choose <b>NM Transaxial</b>	16. Choose <b>NM Transaxial</b>
a. Select an 8 x 8 display grid	a. Select a 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate <b>DELAY TRANS</b> on display screen	c. Annotate <b>DELAY TRANS</b> on display screen
d. <b>SCREENCAP</b> and save as <b>DELAY NM TRANSAXIALS</b>	d. <b>SCREENCAP</b> and save as <b>DELAY NM TRANSAXIALS</b>
17. Repeat step 16 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately	17. Repeat step 16 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately
18. Click <b>File</b> and <b>Save MIP</b>	18. Click <b>File</b> and <b>Save MIP</b>
19. Click <b>File</b> and <b>Save &amp; Exit</b>	19. Click <b>File</b> and <b>Save &amp; Exit</b>
20. Select MIP and rename to <b>LATE MIP</b>	20. Select MIP and rename to <b>LATE MIP</b>
21. Select <b>EARLY</b> and <b>LATE</b> statics	21. Select <b>EARLY</b> and <b>LATE</b> statics
22. Run <b>Parathyroid Imaging</b> processing icon	22. Run <b>Parathyroid Imaging</b> processing icon
23. Adjust intensity accordingly	23. Adjust intensity accordingly
24. <b>SCREENCAP</b> and save as <b>EARLY LATE STATICS</b>	24. <b>SCREENCAP</b> and save as <b>EARLY LATE STATICS</b>
<b>PACS</b> Send all <b>EARLY NM</b> and <b>FUSED SCREENCAPS</b> , all <b>DELAY NM SCREENCAPS</b> , <b>EARLY AND LATE MIPS</b> , <b>CT CORRECTED</b> , <b>EARLY</b> and <b>LATE STATIC SCREENCAP</b> , <b>IRAC OSEM</b> and <b>IROSEM</b> to <b>PACS</b> .	<b>PACS</b> Send all <b>EARLY NM</b> and <b>FUSED SCREENCAPS</b> , all <b>DELAY NM SCREENCAPS</b> , <b>EARLY AND LATE MIPS</b> , <b>CT CORRECTED</b> , <b>EARLY</b> and <b>LATE STATIC SCREENCAP</b> , <b>TOMOEARLY_IRACRR Transaxials</b> and <b>TOMODELAY_IRNC_Transaxials</b> <b>PACS</b> .

**Recon / Reformat Processing and Screen Captures for SPECT ONLY (no CT)**


Infinia 1 with Xeleris 2.x (Rm E)	Infinia 2/3 & NM640 with Xeleris 3.x (RM A, C and TAC)
1. Select patient and the following files	1. Select patient and the following files
a. TOMOearlyHWKY	a. TOMOearlyHWKY
2. Run <b>VOLUMETRIX FOR HAWKEYE PARATHYROID</b> processing icon	2. Run <b>Parathyroid MI</b> processing icon
	a. Click Original
	b. Click Proceed
3. Select <b>NM Transaxials</b>	3. Select <b>NM Transaxials</b>
a. Create an 8 x 8 display grid	a. Create an 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate <b>EARLY TRANS</b> on display screen	c. Annotate <b>EARLY TRANS</b> on display screen
d. <b>SCREENCAP</b> and save as <b>EARLY NM TRANSAXIALS</b>	d. <b>SCREENCAP</b> and save as <b>EARLY NM TRANSAXIALS</b>
4. Repeat step 3 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately	4. Repeat step 3 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately
5. Click <b>File</b> and <b>Save MIP</b>	5. Click <b>File</b> and <b>Save MIP</b>
6. Click <b>File</b> then <b>Save &amp; Exit</b>	6. Click <b>File</b> then <b>Save &amp; Exit</b>
7. Select MIP and rename to <b>EARLY MIP</b>	7. Select MIP and rename to <b>EARLY MIP</b>
8. Select <b>Tomolate</b>	8. Select <b>Tomolate</b>
9. Run <b>VOLUMETRIX FOR HAWKEYE PARATHYROID</b> processing icon	9. Run <b>VOLUMETRIX FOR HAWKEYE PARATHYROID</b> processing icon
10. Click <b>Resume</b>	10. Click <b>Resume</b>
11. Choose <b>NM Transaxial</b>	11. Choose <b>NM Transaxial</b>
a. Select an 8 x 8 display grid	a. Select a 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate <b>DELAY TRANS</b> on display screen	c. Annotate <b>DELAY TRANS</b> on display screen
d. <b>SCREENCAP</b> and save as <b>DELAY NM TRANSAXIALS</b>	d. <b>SCREENCAP</b> and save as <b>DELAY NM TRANSAXIALS</b>
12. Repeat step 11 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately	12. Repeat step 11 by selecting <b>NM CORONALS</b> and <b>NM SAGITTALS</b> . Annotate appropriately
13. Click <b>File</b> and <b>Save MIP</b>	13. Click <b>File</b> and <b>Save MIP</b>
14. Click <b>File</b> and <b>Save &amp; Exit</b>	14. Click <b>File</b> and <b>Save &amp; Exit</b>
15. Select MIP and rename to <b>LATE MIP</b>	15. Select MIP and rename to <b>LATE MIP</b>
16. Select <b>EARLY</b> and <b>LATE</b> statics	16. Select <b>EARLY</b> and <b>LATE</b> statics
17. Run <b>Parathyroid Imaging</b> processing icon	17. Run <b>Parathyroid Imaging</b> processing icon
18. Adjust intensity accordingly	18. Adjust intensity accordingly
19. <b>SCREENCAP</b> and save as <b>EARLY LATE STATICS</b>	19. <b>SCREENCAP</b> and save as <b>EARLY LATE STATICS</b>
<p><b>PACS</b></p> <p>Send all EARLY NM and FUSED SCREENCAPS, all DELAY NM SCREENCAPS, EARLY AND LATE MIPs, CT CORRECTED, EARLY and LATE STATIC SCREENCAP, IRAC OSEM and IROSEM to PACS.</p>	<p><b>PACS</b></p> <p>Send all EARLY NM and FUSED SCREENCAPS, all DELAY NM SCREENCAPS, EARLY AND LATE MIPs, CT CORRECTED, EARLY and LATE STATIC SCREENCAP, TOMOEARLY_IRACRR Transaxials and TOMODELAY_IRNC_Transaxials PACS.</p>

# Screen Cap Samples


## Early Transaxial




## Early Fused Transaxial



## Early Coronal



## Early Sagittal



---

Reviewed By:

Mike Decorah and Dani Brunk

---

Scott B. Perlman, MD, MS  
Chief, Nuclear Medicine

---

Derek Fuerbringer, CNMT  
Manager, Nuclear Medicine  
University Hospital

---

Kandace Nowakowski  
Manager, Nuclear Medicine  
The American Center

---

John Vetter, PhD, DABR  
Medical Physicist

---

Scott Knishka, RPh, BCNP  
Radiopharmacist

---