

Division of Nuclear Medicine Procedure / Protocol University Hospital and The American Center

PARATHYROID SPECT PROCESSING Tc99m
UPDATED: December 2019

CPT CODE: N/A

Recon / Reformat Processing and Screen Captures for SPECT/CT

Reconstruction Options:	Infinia 1 (Rm E)	Infinia 2/3 & NM640 with Xeleris 3.x (RM A, C and TAC)
Reconstruction Type	OSEM/MLEM	OSEM
Number of OSEM Iterations	OSEM 2 MLEM 0	2
Max Number of Subsets	10	10
Correction Options:		
Corrections		Resolutions recovery correction & Attenuation Correction
Additional Set		Resolution Recovery Correction and Attenuation Correction
Filters:		
Pre-Filter	Butterworth	
Pre Critical Frequency	0.50	
Pre Power	10	
Post-Filter	Hann	Butterworth
Critical Frequency	1.2	0.50
Power	10.0	10.0

Infinia 1 with Xeleris 2.x (Rm E)	Infinia 2/3 & NM640 with Xeleris 3.x (RM A, C and TAC)
1. Select the Patient and the following files	1. Select the Patient and the following files
a. TOMOearlyHWKY	a. TOMOearlyHWKY
b. CT TOMOearlyHWKY	b. CT TOMOearlyHWKY
c. ATT MAP TOMOearlyHWKY	c. ATT MAP TOMOearlyHWKY
2. Run VOLUMETRIX FOR HAWKEYE PARATHYROID processing icon	2. Run Parathyroid MI processing icon
	a. Click Original
	b. Click Proceed
	c. Click Passed
3. Select NM Transaxials	3. Select NM Transaxials
	a. Change SUM:1 and STEP:1
a. Create an 8 x 8 display grid	b. Create an 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	c. Center area of interest within the grid and adjust intensity accordingly
c. Annotate EARLY TRANS on display screen	d. Annotate EARLY TRANS on display screen
d. SCREENCAP and save as EARLY NM TRANSAXIALS	e. SCREENCAP and save as EARLY NM TRANSAXIALS

4. Repeat step 3 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately	4. Repeat step 3 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately
5. Select FUSED Transaxials	5. Select FUSED Transaxials
a. Create an 8 x 8 display grid	a. Create an 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate EARLY TRANS on display screen	c. Annotate EARLY TRANS on display screen
d. SCREENCAP and save as EARLY FUSED TRANSAXIALS	d. SCREENCAP and save as EARLY FUSED TRANSAXIALS
6. Repeat step 5 by selecting FUSED CORONALS AND FUSED SAGITTALS . Annotate appropriately	6. Repeat step 5 by selecting FUSED CORONALS AND FUSED SAGITTALS . Annotate appropriately
7. Click File and Save MIP	7. Click File and Save MIP
8. Click File then Save & Exit	8. Click File then Save & Exit
9. Select MIP and rename to EARLY MIP	9. Select MIP and rename to EARLY MIP
10. Select CT TOMOearlyHWKY	10. Select CT TOMOearlyHWKY
11. Run Convert CT to Hounsfield Units processing icon	11. Run Convert CT to Hounsfield Units processing icon
12. Click File and Exit	12. Click File and Exit
13. Select Tomolate	13. Select Tomolate
14. Run VOLUMETRIX FOR HAWKEYE PARATHYROID processing icon	14. Run VOLUMETRIX FOR HAWKEYE PARATHYROID processing icon
15. Click Resume	15. Click Resume
16. Choose NM Transaxial	16. Choose NM Transaxial
a. Select an 8 x 8 display grid	a. Select a 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate DELAY TRANS on display screen	c. Annotate DELAY TRANS on display screen
d. SCREENCAP and save as DELAY NM TRANSAXIALS	d. SCREENCAP and save as DELAY NM TRANSAXIALS
17. Repeat step 16 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately	17. Repeat step 16 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately
18. Click File and Save MIP	18. Click File and Save MIP
19. Click File and Save & Exit	19. Click File and Save & Exit
20. Select MIP and rename to LATE MIP	20. Select MIP and rename to LATE MIP
21. Select EARLY and LATE statics	21. Select EARLY and LATE statics
22. Run Parathyroid Imaging processing icon	22. Run Parathyroid Imaging processing icon
23. Adjust intensity accordingly	23. Adjust intensity accordingly
24. SCREENCAP and save as EARLY LATE STATICS	24. SCREENCAP and save as EARLY LATE STATICS
PACS Send all EARLY NM and FUSED SCREENCAPS , all DELAY NM SCREENCAPS , EARLY AND LATE MIPS , CT CORRECTED , EARLY and LATE STATIC SCREENCAP , IRAC OSEM and IROSEM to PACS .	PACS Send all EARLY NM and FUSED SCREENCAPS , all DELAY NM SCREENCAPS , EARLY AND LATE MIPS , CT CORRECTED , EARLY and LATE STATIC SCREENCAP , TOMOEARLY_IRACRR Transaxials and TOMODELAY_IRNC_Transaxials PACS .

Recon / Reformat Processing and Screen Captures for SPECT ONLY (no CT)


Infinia 1 with Xeleris 2.x (Rm E)	Infinia 2/3 & NM640 with Xeleris 3.x (RM A, C and TAC)
1. Select patient and the following files	1. Select patient and the following files
a. TOMOearlyHWKY	a. TOMOearlyHWKY
2. Run VOLUMETRIX FOR HAWKEYE PARATHYROID processing icon	2. Run Parathyroid MI processing icon
	a. Click Original
	b. Click Proceed
3. Select NM Transaxials	3. Select NM Transaxials
a. Create an 8 x 8 display grid	a. Create an 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate EARLY TRANS on display screen	c. Annotate EARLY TRANS on display screen
d. SCREENCAP and save as EARLY NM TRANSAXIALS	d. SCREENCAP and save as EARLY NM TRANSAXIALS
4. Repeat step 3 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately	4. Repeat step 3 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately
5. Click File and Save MIP	5. Click File and Save MIP
6. Click File then Save & Exit	6. Click File then Save & Exit
7. Select MIP and rename to EARLY MIP	7. Select MIP and rename to EARLY MIP
8. Select Tomolate	8. Select Tomolate
9. Run VOLUMETRIX FOR HAWKEYE PARATHYROID processing icon	9. Run VOLUMETRIX FOR HAWKEYE PARATHYROID processing icon
10. Click Resume	10. Click Resume
11. Choose NM Transaxial	11. Choose NM Transaxial
a. Select an 8 x 8 display grid	a. Select a 8 x 8 display grid
b. Center area of interest within the grid and adjust intensity accordingly	b. Center area of interest within the grid and adjust intensity accordingly
c. Annotate DELAY TRANS on display screen	c. Annotate DELAY TRANS on display screen
d. SCREENCAP and save as DELAY NM TRANSAXIALS	d. SCREENCAP and save as DELAY NM TRANSAXIALS
12. Repeat step 11 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately	12. Repeat step 11 by selecting NM CORONALS and NM SAGITTALS . Annotate appropriately
13. Click File and Save MIP	13. Click File and Save MIP
14. Click File and Save & Exit	14. Click File and Save & Exit
15. Select MIP and rename to LATE MIP	15. Select MIP and rename to LATE MIP
16. Select EARLY and LATE statics	16. Select EARLY and LATE statics
17. Run Parathyroid Imaging processing icon	17. Run Parathyroid Imaging processing icon
18. Adjust intensity accordingly	18. Adjust intensity accordingly
19. SCREENCAP and save as EARLY LATE STATICS	19. SCREENCAP and save as EARLY LATE STATICS
<p>PACS</p> <p>Send all EARLY NM and FUSED SCREENCAPS, all DELAY NM SCREENCAPS, EARLY AND LATE MIPs, CT CORRECTED, EARLY and LATE STATIC SCREENCAP, IRAC OSEM and IROSEM to PACS.</p>	<p>PACS</p> <p>Send all EARLY NM and FUSED SCREENCAPS, all DELAY NM SCREENCAPS, EARLY AND LATE MIPs, CT CORRECTED, EARLY and LATE STATIC SCREENCAP, TOMOEARLY_IRACRR Transaxials and TOMODELAY_IRNC_Transaxials PACS.</p>

Screen Cap Samples


Early Transaxial




Early Fused Transaxial




Early Coronal




Early Fused Coronal




Early Sagittal



Early Fused Sagittal



Statics Early and Delayed



Reviewed By:

Mike Decorah and Dani Brunk

Scott B. Perlman, MD, MS
Chief, Nuclear Medicine

Derek Fuerbringer, CNMT
Manager, Nuclear Medicine
University Hospital

Kandace Nowakowski
Manager, Nuclear Medicine
The American Center

John Vetter, PhD, DABR
Medical Physicist

Scott Knishka, RPh, BCNP
Radiopharmacist
