

Shoulder TUMOR (Adult) SIGNA Artist

Technical Instructions

--power inj/subtractions--

Lt 3PI Loc

Rt 3PI Loc

Ax T2 Fat DL

Ax T1 DL

Cor T1 DL

Cor T2 Fat DL

Sag T2 Fat DL

Ax T1 Lava-Flex Pre

Ax T1 Lava-Flex 2min

Cor T1 Lava-Flex

Sag T1 Lava-Flex

+c Ax T1 Fat DL

Lt 3PI Loc

Details tab

Protocol notes

Lt 3PI Loc

PATIENT POSITION

Coil Configuration:	Inv16chShoulder
Plane:	3-PLANE

IMAGING PARAMETERS

Imaging Mode:	2D
Imaging Options:	Seq, EDR, TRF, Fast, SS, ARC
Pulse Sequence:	Spin Echo
ARC Phase:	2.00

Image Recon

Recon DL Strength:	Off
Intensity Correction:	NONE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	32.0
Slice Thickness:	6.00
Slice Spacing:	0.0
Overlap Locations:	0
No Phase Wrap:	1.50

ACQ TIMING

Acq Voxel Size:	1.3 X 2.0 X 6.0
Freq:	256
Phase:	160
Freq DIR:	Unswap
# of Acq. Before Pause:	0
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective

SCAN TIMING

TE:	100.0
Number of Echoes:	1
TR:	Minimum
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	No
------------------	----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

USER CVS

User CV1:	1.00
-----------	------

Rt 3PI Loc

Details tab

Protocol notes

Rt 3PI Loc

PATIENT POSITION

Coil Configuration:	Inv16chShoulder
Plane:	3-PLANE

IMAGING PARAMETERS

Imaging Mode:	2D
Imaging Options:	Seq, EDR, TRF, Fast, SS, ARC
Pulse Sequence:	Spin Echo
ARC Phase:	2.00

Image Recon

Recon DL Strength:	Off
Intensity Correction:	NONE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	32.0
Slice Thickness:	6.00
Slice Spacing:	5.0
Overlap Locations:	0
No Phase Wrap:	1.50

ACQ TIMING

Acq Voxel Size:	1.3 X 2.0 X 6.0
Freq:	256
Phase:	160
Freq DIR:	Unswap
# of Acq. Before Pause:	0
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective

SCAN TIMING

TE:	100.0
Number of Echoes:	1
TR:	Minimum
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	No
------------------	----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

USER CVS

User CV1:	1.00
-----------	------

Ax T2 Fat DL

Details tab

Protocol notes

Ax T2 Fat DL

PATIENT POSITION

Coil Configuration: GE40chGEMPA, Inv16chCardiacAA
Plane: OBLIQUE

IMAGING PARAMETERS

Imaging Mode: 2D
Imaging Options: EDR, TRF, Fast
Pulse Sequence: FSE-XL

Image Recon

Recon DL Strength: Medium
Intensity Correction: PURE
Intensity Filter: None
Save Original: Off

SCAN RANGE

FOV: 16.0
Slice Thickness: 5.00
Slice Spacing: 1.5
Overlap Locations: 0
No Phase Wrap: 1.50

ACQ TIMING

Acq Voxel Size:	0.6 X 0.6 X 5.0
Freq:	288
Phase:	256
Freq DIR:	Unswap
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	On
RF Drive Mode:	Quadrature
Excitation Mode:	Selective
Fat Shift DIR:	Normal

SCAN TIMING

Flip Angle:	160.0
TE:	60.0
Number of Echoes:	1
TR:	2149.0
Echo Train Length:	8
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	No
------------------	----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

SAT

Fat/Water Saturation:

Fat

USER CVS

User CV6:	None
User CV9:	None
User CV13:	None
User CV15:	None
User CV16:	None
User CV18:	None
User CV19:	1.00
User CV21:	None
User CV26:	1.00
TR Min:	3000.0
TR Max:	10000.0

Non-Viewable Parameters

ñ

Ax T1 DL

Details tab

Protocol notes

Ax T1 DL

PATIENT POSITION

Coil Configuration:	GE40chGEMPA, Inv16chCardiacAA
Plane:	OBLIQUE

IMAGING PARAMETERS

Imaging Mode:	2D
Imaging Options:	EDR, TRF, Fast, ARC
Pulse Sequence:	FSE-XL
ARC Phase:	1.50

Image Recon

Recon DL Strength:	Medium
Intensity Correction:	SCENIC
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	16.0
Slice Thickness:	5.00
Slice Spacing:	1.5
Overlap Locations:	0
No Phase Wrap:	1.50

ACQ TIMING

Acq Voxel Size:	0.5 X 0.6 X 5.0
Freq:	320
Phase:	288
Freq DIR:	Unswap
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective
Fat Shift DIR:	Normal

SCAN TIMING

Flip Angle:	160.0
TE:	20.0
Number of Echoes:	1
TR:	478.0
Echo Train Length:	4
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	No
------------------	----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

User CV6:	None
User CV9:	None
User CV15:	None
User CV16:	None
User CV18:	None
User CV21:	None
User CV26:	1.00
TR Min:	300.0
TR Max:	750.0

Non-Viewable Parameters

ñ

Cor T1 DL

Details tab

Protocol notes

Cor T1 DL

PATIENT POSITION

Coil Configuration:	GE40chGEMPA, Inv16chCardiacAA
Plane:	OBLIQUE

IMAGING PARAMETERS

Imaging Mode:	2D
Imaging Options:	EDR, TRF, Fast, ARC
Pulse Sequence:	FSE-XL
ARC Phase:	1.50

Image Recon

Recon DL Strength:	Medium
Intensity Correction:	SCENIC
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	18.0
Slice Thickness:	4.00
Slice Spacing:	1.0
Overlap Locations:	0
No Phase Wrap:	1.50

ACQ TIMING

Acq Voxel Size:	0.6 X 0.6 X 4.0
Freq:	320
Phase:	288
Freq DIR:	Unswap
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	On
RF Drive Mode:	Quadrature
Excitation Mode:	Selective
Fat Shift DIR:	Normal

SCAN TIMING

Flip Angle:	160.0
TE:	20.0
Number of Echoes:	1
TR:	511.0
Echo Train Length:	4
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	No
------------------	----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

User CV6:	None
User CV9:	None
User CV15:	None
User CV16:	None
User CV18:	None
User CV21:	None
User CV26:	1.00

Non-Viewable Parameters

ñ

Cor T2 Fat DL

Details tab

Protocol notes

Cor T2 Fat DL

PATIENT POSITION

Coil Configuration:	GE40chGEMPA, Inv16chCardiacAA
Plane:	OBLIQUE

IMAGING PARAMETERS

Imaging Mode:	2D
Imaging Options:	EDR, TRF, Fast
Pulse Sequence:	FSE-XL

Image Recon

Recon DL Strength:	Medium
Intensity Correction:	PURE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	18.0
Slice Thickness:	4.00
Slice Spacing:	1.0
Overlap Locations:	0
No Phase Wrap:	1.50

ACQ TIMING

Acq Voxel Size:	0.6 X 0.7 X 4.0
Freq:	288
Phase:	256
Freq DIR:	Unswap
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	On
RF Drive Mode:	Quadrature
Excitation Mode:	Selective
Fat Shift DIR:	Normal

SCAN TIMING

Flip Angle:	160.0
TE:	60.0
Number of Echoes:	1
TR:	3000.0
Echo Train Length:	8
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	No
------------------	----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

SAT

Fat/Water Saturation:

Fat

USER CVS

User CV6:	None
User CV9:	None
User CV13:	None
User CV15:	None
User CV16:	None
User CV18:	None
User CV19:	1.00
User CV21:	None
User CV26:	1.00

Non-Viewable Parameters

ñ

Sag T2 Fat DL

Details tab

Protocol notes

Sag T2 Fat DL

PATIENT POSITION

Coil Configuration:	GE40chGEMPA, Inv16chCardiacAA
Plane:	OBLIQUE

IMAGING PARAMETERS

Imaging Mode:	2D
Imaging Options:	EDR, TRF, Fast
Pulse Sequence:	FSE-XL

Image Recon

Recon DL Strength:	Medium
Intensity Correction:	PURE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	18.0
Slice Thickness:	4.00
Slice Spacing:	1.0
Overlap Locations:	0
No Phase Wrap:	1.50

ACQ TIMING

Acq Voxel Size:	0.6 X 0.7 X 4.0
Freq:	288
Phase:	256
Freq DIR:	Unswap
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	On
RF Drive Mode:	Quadrature
Excitation Mode:	Selective
Fat Shift DIR:	Normal

SCAN TIMING

Flip Angle:	160.0
TE:	60.0
Number of Echoes:	1
TR:	3000.0
Echo Train Length:	8
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	No
------------------	----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

SAT

Fat/Water Saturation:

Fat

USER CVS

User CV6:	None
User CV9:	None
User CV13:	None
User CV15:	None
User CV16:	None
User CV18:	None
User CV19:	1.00
User CV21:	None
User CV26:	1.00

Non-Viewable Parameters

ñ

Ax T1 Lava-Flex Pre

Details tab

Protocol notes

Ax T1 Lava-Flex Pre

PATIENT POSITION

Coil Configuration:	GE30chAirAA, GE40chGEMPA
Plane:	AXIAL

IMAGING PARAMETERS

Imaging Mode:	3D
Imaging Options:	EDR, Fast, ARC, Flex
Pulse Sequence:	LAVA
IDEAL:	15
ARC Phase:	2.00
ARC Slice:	1.00

Image Recon

Recon DL Strength:	Off
Intensity Correction:	PURE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	26.0
Slice Thickness:	1.50
Slices per Slab:	172
Overlap Locations:	0

No Slab Wrap:	1.04
No Phase Wrap:	1.00
Slice Resolution(%):	50.00

ACQ TIMING

Acq Voxel Size:	1.0 X 1.2 X 3.0
Freq:	260
Phase:	220
Freq DIR:	R/L
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective

SCAN TIMING

Flip Angle:	15.0
TE:	Min Full
Number of Echoes:	2
Receiver Bandwidth:	90.91

CONTRAST

Contrast Yes/No:	Yes
------------------	-----

OTHERS

3D Geometry Correction:	On
-------------------------	----

Anatomical Region:

Rectum

USER CVS

User CV4:	None
User CV6:	1.00
User CV13:	None
User CV16:	None
User CV30:	1.00
User CV34:	70.00

Non-Viewable Parameters

ñ

Ax T1 Lava-Flex 2min

Details tab

Protocol notes

Ax T1 Lava-Flex 2min

PATIENT POSITION

Coil Configuration:	GE30chAirAA, GE40chGEMPA
Plane:	AXIAL

IMAGING PARAMETERS

Imaging Mode:	3D
Imaging Options:	EDR, Fast, MPh, ARC, Flex
Pulse Sequence:	LAVA
IDEAL:	3
ARC Phase:	2.00
ARC Slice:	1.00

Image Recon

Recon DL Strength:	Off
Intensity Correction:	PURE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	26.0
Slice Thickness:	1.50
Slices per Slab:	172
Overlap Locations:	0

No Slab Wrap:	1.04
No Phase Wrap:	1.00
Slice Resolution(%):	50.00

ACQ TIMING

Acq Voxel Size:	1.0 X 1.2 X 3.0
Freq:	260
Phase:	220
Freq DIR:	R/L
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective

SCAN TIMING

Flip Angle:	15.0
TE:	Min Full
Number of Echoes:	2
Receiver Bandwidth:	90.91

MULTI-PHASE

Slice per Location:	1
Seperate Series:	0
Trigger Delay without AV:	0
Mask Phase:	Off

Mask Pause: **Off**

CONTRAST

Contrast Yes/No: Yes

OTHERS

3D Geometry Correction: On

Anatomical Region: Rectum

USER CVS

User CV4: None

User CV6: 1.00

User CV13: None

User CV16: None

User CV30: 1.00

User CV34: 70.00

Non-Viewable Parameters

ñ

Cor T1 Lava-Flex

Details tab

Protocol notes

Cor T1 Lava-Flex

PATIENT POSITION

Coil Configuration:	GE30chAirAA, GE40chGEMPA
Plane:	CORONAL

IMAGING PARAMETERS

Imaging Mode:	3D
Imaging Options:	EDR, Fast, ARC, Flex
Pulse Sequence:	LAVA
IDEAL:	3
ARC Phase:	2.00
ARC Slice:	1.00

Image Recon

Recon DL Strength:	Off
Intensity Correction:	PURE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	26.0
Slice Thickness:	1.50
Slices per Slab:	172
Overlap Locations:	0

No Slab Wrap:	1.04
No Phase Wrap:	1.00
Slice Resolution(%):	50.00

ACQ TIMING

Acq Voxel Size:	1.0 X 1.2 X 3.0
Freq:	260
Phase:	220
Freq DIR:	S/I
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective

SCAN TIMING

Flip Angle:	15.0
TE:	Min Full
Number of Echoes:	2
Receiver Bandwidth:	90.91

CONTRAST

Contrast Yes/No:	Yes
------------------	-----

OTHERS

3D Geometry Correction:	On
-------------------------	----

Anatomical Region:

Rectum

USER CVS

User CV4:	None
User CV6:	1.00
User CV13:	None
User CV16:	None
User CV30:	1.00
User CV34:	70.00

Non-Viewable Parameters

ñ

Sag T1 Lava-Flex

Details tab

Protocol notes

Sag T1 Lava-Flex

PATIENT POSITION

Coil Configuration:	GE30chAirAA, GE40chGEMPA
Plane:	SAGITTAL

IMAGING PARAMETERS

Imaging Mode:	3D
Imaging Options:	EDR, Fast, ARC, Flex
Pulse Sequence:	LAVA
IDEAL:	3
ARC Phase:	2.00
ARC Slice:	1.00

Image Recon

Recon DL Strength:	Off
Intensity Correction:	PURE
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	26.0
Slice Thickness:	1.50
Slices per Slab:	172
Overlap Locations:	0

No Slab Wrap:	1.04
No Phase Wrap:	1.00
Slice Resolution(%):	50.00

ACQ TIMING

Acq Voxel Size:	1.0 X 1.2 X 3.0
Freq:	260
Phase:	220
Freq DIR:	S/I
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective

SCAN TIMING

Flip Angle:	15.0
TE:	Min Full
Number of Echoes:	2
Receiver Bandwidth:	90.91

CONTRAST

Contrast Yes/No:	Yes
------------------	-----

OTHERS

3D Geometry Correction:	On
-------------------------	----

Anatomical Region:

Rectum

USER CVS

User CV4:	None
User CV6:	1.00
User CV13:	None
User CV16:	None
User CV30:	1.00
User CV34:	70.00

Non-Viewable Parameters

ñ

+c Ax T1 Fat DL

Details tab

Protocol notes

+c Ax T1 Fat DL

PATIENT POSITION

Coil Configuration:	GE40chGEMPA, Inv16chCardiacAA
Plane:	OBLIQUE

IMAGING PARAMETERS

Imaging Mode:	2D
Imaging Options:	EDR, TRF, Fast, ARC
Pulse Sequence:	FSE-XL
ARC Phase:	1.50

Image Recon

Recon DL Strength:	Medium
Intensity Correction:	SCENIC
Intensity Filter:	None
Save Original:	Off

SCAN RANGE

FOV:	16.0
Slice Thickness:	5.00
Slice Spacing:	1.5
Overlap Locations:	0
No Phase Wrap:	1.50

ACQ TIMING

Acq Voxel Size:	0.5 X 0.6 X 5.0
Freq:	320
Phase:	288
Freq DIR:	Unswap
NEX:	1.00
Phase FOV:	1.00
Auto Shim:	Auto
Phase Correct:	Off
RF Drive Mode:	Quadrature
Excitation Mode:	Selective
Fat Shift DIR:	Normal

SCAN TIMING

Flip Angle:	160.0
TE:	20.0
Number of Echoes:	1
TR:	478.0
Echo Train Length:	4
Receiver Bandwidth:	62.50

CONTRAST

Contrast Yes/No:	Yes
------------------	-----

OTHERS

Anatomical Region:	Shoulder
--------------------	----------

Fat/Water Saturation:

Fat

USER CVS

User CV6:	None
User CV9:	None
User CV13:	None
User CV15:	None
User CV16:	None
User CV18:	None
User CV19:	1.00
User CV21:	None
User CV26:	1.00
TR Min:	300.0
TR Max:	750.0

Non-Viewable Parameters

ñ