CALF Proximal PERONEAL NERVE NEUROGRAM 11-15-16 Monitored	
 *Call rad in A.M. to have them review case. That way coil selection and pathology can be reviewed 3T ONLY 8 Channel Cardiac Coil 1. 3 Pl loc 2. Ax T1 3. Ax T2 dk fat (Separate into upper and lower half) 4. Sag T1 5. Sag T2 IDEAL 6. Coronal T1 7. Cor T2 dk fat Give .1mmol/kg Multihance Low eGFR Inpatient Dose: No Change: 3 planes +c T1 dk fat 	
Nerve	<image/> <caption><caption><image/></caption></caption>