

Division of Nuclear Medicine Procedure / Protocol

RADIOPHARMACY REFRIGERATOR (CLIA-88)
UPDATED: AUGUST 2011

CPT CODE: N/A

Purpose: In order to ensure storage of radiopharmaceuticals at manufacturers' recommended temperature range.

Policy: The temperature check will be done daily, except for Saturdays, Sundays, and Holidays.

<u>Temperature Range</u>	<u>Item</u>
Freezer	< 0° C
Refrigerator	2-8° C

Action Plan: If temperature is out of range, it will be checked again later in the day. If a problem persists, Plant Engineering will be notified.

WEEKLY REVIEW.

Reviewed By: S. Perlman, D. Fuerbringer

Scott B. Perlman, MD, MS
Chief, Nuclear Medicine

Derek Fuerbringer, CNMT
Manager, Nuclear Medicine